HIV-1/2 Ag/Ab Combo

Key to symbols used / Schlüssel für die verwendeten Symbole / Clave de los símbolos utilizados / Légende des symboles utilisés / Legenda dei simboli utilizzati / Descrição dos símbolos utilizados

REF Catalogue Number / Katalognummer / Número de catálogo / Référence catalogue / Numero di catalogo / N° de Catálogo

Contains Sufficient for 100 tests / Enthält eine ausreichende Menge für 100 Tests / Contiene material suficiente para realizar 100 pruebas / Permet de réaliser 100 tests / Contiene reagenti sufficienti per 100 tests / Contém o suficiente para 100 testes

IVD In Vitro Diagnostic Medical Device / In Vitro - Diagnostika / Dispositivo de diagnóstico médico *in vitro* / Dispositif médical de diagnostic *In Vitro* / Dispositivo medico diagnostico *In Vitro* / Dispositivo Médico para Diagnóstico *In Vitro*

Store at 2-30°C / Bei 2-30°C lagern / Guardar a temperaturas entre 2 y 30°C / Conserver entre 2-30°C / Conservare a 2-30°C / Armazenar a 2-30°C

**Determine™ HIV-1/2 Ag/Ab Combo**

EN

This package insert must be read carefully prior to use. Package insert instructions must be followed accordingly. Reliability of assay results cannot be guaranteed if there are deviations from the instructions in this package insert.

NAME AND INTENDED USE

Alere Determine™ HIV-1/2 Ag/Ab Combo is an *In Vitro*, visually read, qualitative immunoassay for the simultaneous detection of free non immunocomplexed HIV-1 p24 antigen (Ag) and antibodies (Ab) to HIV-1 and HIV-2 in human blood. The test specimen can be serum, plasma, fingerstick or venous whole blood. The test is intended as an aid to detect HIV-1 p24 antigen and antibodies to HIV-1/HIV-2 from infected individuals.

SUMMARY AND EXPLANATION OF THE TEST

AIDS (Acquired Immunodeficiency Syndrome) is characterized by changes in the population of T-cell lymphocytes. In an infected individual, the virus causes depletion of helper T-cells, which leaves the person susceptible to opportunistic infections and some malignancies.

The virus that causes AIDS exists as two related types known as HIV-1 and HIV-2. The presence of HIV first elicits the secretion of p24 antigen^{1,2} followed by the production of specific antibodies to either HIV-1 or HIV-2.^{3,4,5}

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

Alere Determine™ HIV-1/2 Ag/Ab Combo is an immunochromatographic test for the qualitative detection of p24 antigen and antibodies to HIV-1 and HIV-2.

Specimen is added to the sample pad. The specimen mixes with a biotinylated anti-p24 antibody, selenium colloid-antigen conjugate and selenium colloid-anti-p24 antibody. This mixture continues to migrate through the solid phase to the immobilized avidin, recombinant antigens and synthetic peptides at the patient window sites.

If antibodies to HIV-1 and/or HIV-2 are present in the specimen, the antibodies bind to the antigen-selenium colloid and to the immobilized recombinant antigens and synthetic peptides, forming one red bar at the patient HIV Antibody window site. If antibodies to HIV-1 and/or HIV-2 are absent the antigen-selenium colloid flows past the patient window, and no red bar is formed at the patient HIV Antibody window site.

If free non immunocomplexed HIV-1 p24 antigen (Ag), is present in the specimen, the antigen binds to the biotinylated anti-p24 from the sample pad and the selenium colloid anti-p24 antibody and it binds to an immobilized avidin forming a red bar at the patient HIV Antigen window site. If p24 antigen is not present both the biotinylated anti-p24 and selenium colloid anti-p24 antibody flow past the patient window, and no red bar is formed at the patient HIV Antigen window site.

To ensure assay validity, a procedural control bar is incorporated in the assay device.

CONTENTS

Alere Determine™ HIV-1/2 Ag/Ab Combo 100 Test (7D2643)

- **Alere Determine™ HIV-1/2 Ag/Ab Combo** Test Card, 10 cards (10 tests/card) coated with HIV-1/2 recombinant antigen and synthetic peptides, anti-p24 antibodies and avidin.

ACCESSORIES (required but not provided)

For testing Whole Blood samples

1 Bottle (2.5 mL) Chase Buffer (7D2243) prepared in phosphate buffer. Preservatives: Antimicrobial Agents.

Whole Blood (fingerstick assay)

EDTA Capillary Tube (7D2222)

Lancet

WARNINGS AND PRECAUTIONS

For *In Vitro* Diagnostic Use.

For Professional use only.

Safety data sheet available for professional user on request.

Patients with elevated triglycerides levels may test false reactive with the **Alere Determine™ HIV-1/2 Ag/Ab Combo**.

CAUTION:

Appropriate biosafety practices^{7,8} should be used when handling specimens and reagents. These precautions include, but are not limited to the following:

- Wear gloves.
- Do not pipette by mouth.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in areas where these materials are handled.
- Clean and disinfect all spills of specimens or reagents using suitable disinfectant, such as 0.5% sodium hypochlorite, or other suitable disinfectant.⁶
- Decontaminate and dispose of all specimens, reagents, and other potentially contaminated materials in accordance with local regulations.^{7,8}

STORAGE

Alere Determine™ HIV-1/2 Ag/Ab Combo Test Cards and Chase Buffer must be stored at 2-30°C until expiration date.

- Kit components are stable until expiration date when handled and stored as directed. Do not use kit components beyond expiration date.
- Immediately reseal all unused tests in the foil pouch containing the desiccant by pressing seal from end to end to close.
- Do not use devices that have become wet or if the packaging has become damaged.

SPECIMEN COLLECTION

Serum, Plasma, and Whole Blood Collection by Venipuncture.

Human serum, plasma, and whole blood collected by venipuncture should be collected aseptically in such a way as to avoid hemolysis.

Separate the serum from the clot or plasma from the packed cells as soon as possible to avoid any hemolysis.

NOTE: For whole blood and plasma specimens, EDTA collection tubes must be used.

Whole Blood Collection by Fingerstick⁹

Before collecting a fingerstick specimen, place an EDTA capillary tube on a clean dry surface.

1. Choose the fingertip of the middle, ring, or index finger (whichever is the least calloused). Warm the hand as needed with a warm, moist towel or warm water to increase blood flow.
2. Clean fingertip with alcohol; allow to air dry.
3. Position the hand palm-side up. Place the lancet off-center on the fingertip. Firmly press the lancet against the finger and puncture the skin. Dispose of the lancet in an appropriate biohazard sharps container.
4. Wipe away the first drop of blood with a sterile gauze pad.
5. Hold the finger lower than the elbow and apply gentle, intermittent pressure to the base of the punctured finger several times. Touch the tip of the EDTA Capillary Tube to the drop of blood*. Avoid air bubbles.

*If EDTA Capillary Tubes (No. 7D2222) will be used, fill the tube with blood between the 2 marked lines (50µL).

**SPECIMEN STORAGE**

Serum and plasma specimens should be stored at 2-8°C if the test is to be run within 7 days of collection. If testing is delayed more than 7 days, the specimen should be frozen (-20°C or colder).

Avoid repeated freeze/thaw cycles. Specimens that have been frozen and thawed more than 3 times cannot be used.

All frozen specimens must be centrifuged at 10,000g for 5 min at room temperature. Carefully remove the 50 µL test sample from the supernatant. If a lipid layer is formed on the surface of the liquid, ensure that the sample is taken from the clear liquid below that layer.

Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 7 days of collection. Do not freeze whole blood specimens.

If stored at -20°C, bring specimen to room temperature before testing. Mix specimen well by gentle inversion of the tube immediately before testing.

Whole blood collected by fingerstick should be tested immediately.

TEST PROCEDURE

The desired number of test units from the 10-test card can be removed by bending and tearing at the perforation.

NOTE:

• Removal of the test units should start from the right side of the test card to preserve the lot number which appears on the left side of the test card.

• Assay should be initiated within 2 hours after removing the protective foil cover from each test.

1. Remove the protective foil cover from each test.

2. For serum or plasma samples:

a. Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).

b. Wait a minimum of 20 minutes from addition of the sample (up to 30 minutes maximum) and read result.

3. For whole blood (venipuncture) samples:

a. Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).

b. Wait one minute, then apply one drop of Chase Buffer to the sample pad.

c. Wait a minimum of 20 minutes from addition of the chase buffer (up to 30 minutes maximum) and read result.

4. For whole blood (fingerstick) samples using an EDTA capillary tube:

a. Place the capillary tube containing the blood sample onto the sample pad (marked by the arrow symbol).

b. Wait until all the blood is transferred from the capillary tube to the sample pad and then immediately apply one drop of Chase Buffer to the sample pad. Caution: do not lift the capillary from the sample pad before all the blood has been transferred – a bubble may form which will prevent the complete transfer of sample and invalidate the test.


c. Wait a minimum of 20 minutes from addition of the chase buffer (up to 30 minutes maximum) and read result.

QUALITY CONTROL

To ensure assay validity, a procedural control is incorporated in the device and is labeled "Control". Any visible line (even very faint) in the control window should be interpreted as a valid result. If the control bar does not turn red by assay completion, the test result is invalid and the sample should be retested.

INTERPRETATION OF RESULTS**ANTIBODY POSITIVE (Two Bars - Control and Ab Bars)**

Red bars appear in both the control window (labeled "Control") and in the Ab bar window (labeled "Ab") of the strip. Any visible red (or grey-red) bar in the patient window should be interpreted as positive.

**ANTIGEN (p24) POSITIVE (Two Bars - Control and Ag Bars)**


Red bars appear in the Control window (labeled "Control"), the Ab bar window (labeled "Ab") and Ag bar window (labeled "Ag") of the strip. Any visible red (or grey-red) bar in the Ab bar and Ag bar windows should be interpreted as positive. The presence of an antigen response suggests that the infection is at an early stage. Follow up testing may be suggested in order to track the expected future detection of antibodies.

ANTIBODY POSITIVE AND ANTIGEN (p24) POSITIVE (Three Bars - Control, Ab and Ag Bars)

Red bars appear in the Control window (labeled "Control"), the Ab bar window (labeled "Ab") and Ag bar window (labeled "Ag") of the strip. Any visible red (or grey-red) bar in the Ab bar and Ag bar windows should be interpreted as positive. The presence of an antigen response suggests that the infection is at an early stage.

NEGATIVE (One Bar)

One red bar appears in the control window of the strip (labeled "Control"), and no red bar appears in the patient windows of the strip (labeled "Ag" and "Ab").

**INVALID (No Bar)**

If there is no red bar in the control window of the strip, and even if a red bar appears in one of the patient windows of the strip, the result is invalid and should be repeated.

• The test result is positive even if the patient bars appear lighter or darker than the control bar.

• If an invalid result occurs repeatedly, or for technical assistance, contact your local distributor or call Technical Support.

LIMITATIONS OF THE PROCEDURE

- **Alere Determine™ HIV-1/2 Ag/Ab Combo** is designed to simultaneously detect antibodies to HIV-1 and/or HIV-2 and free non immunocomplexed HIV-1 p24 antigen (Ag), in human serum, plasma and whole blood. Other body fluids or pooled specimens may not give accurate results and should not be used.
- The intensity of the Ab and Ag bars does not correlate to the titer of antibody and antigen in the specimen.
- No test provides absolute assurance that a specimen does not contain low levels of HIV-1 p24 antigen and/or antibodies to HIV-1 and HIV-2 such as those present at a very early stage of infection.
- A negative result for both antibodies to HIV and p24 antigen does not preclude the possibility of exposure to or infection with HIV-1 or HIV-2 viruses.
- The test is designed to increase the HIV case load by improving the detection of HIV infections vs other rapid tests and 3rd generation ELISA. The sensitivity of **Alere Determine™ HIV-1/2 Ag/Ab Combo** is not equivalent to 4th generation HIV ELISA, p24 EIA or PCR, which will limit its ability to detect acute infections in comparison to those methods.
- A false negative result could occur in early infections with some low prevalence non-B HIV subtypes which at that stage are less detectable by the **Alere Determine™ HIV-1/2 Ag/Ab Combo** assay.
- Where clinical presentation or other data would suggest an inconsistent result then the patient should be tested by PCR immediately and/or retested for antibodies to HIV >21 days after the original testing.
- A positive result for antibodies to HIV with a negative result for p24 antigen does not preclude the possibility of acute infection.
- Positive results should be confirmed using another method and the results should be evaluated in light of the overall clinical evaluation before a diagnosis is made.
- No significant interference was found up to EDTA-2Na (200 mg/dL), Hemoglobin (500 mg/dL), Bilirubin (20 mg/dL) and Triglyceride (1 g/dL).
- Infants born to HIV-infected mothers may carry maternal antibodies and will test antibody positive until eighteen months of age, which may not necessarily indicate the true infection status of the new born. The use of p24 antigen testing to exclude infection in neonates (up to around eighteen months) is not recommended by CDC, because of poor sensitivity, especially in the presence of HIV antibody. Definitive diagnosis of HIV infection in early infancy requires other assays, including HIV nucleic acid test or viral culture.⁹
- Some known HIV-infected persons taking Antiretroviral medication have been shown to produce false negative results when tested by rapid kits.^{10,11,12}

PERFORMANCE CHARACTERISTICS

The performance of **Alere Determine™ HIV-1/2 Ag/Ab Combo** has been determined by testing specimens from random blood donors, from patients with HIV infection, patients at risk of HIV infection or in other clinical categories and commercial seroconversion panels. In addition, 30 fresh matched samples confirmed positive for HIV-1 were tested. The performance evaluations were conducted in nine clinical studies in Europe, Africa, Asia and South America.

SENSITIVITY

Sensitivity has been evaluated by testing confirmed HIV Ab positive samples, commercial seroconversion panels and specimens from primary (acute) HIV infected patients.

1. HIV Antibody positive specimens:**Table I**
Sensitivity of Alere Determine™ HIV-1/2 Ag/Ab Combo

Types	Number of Specimens Tested	Positive by Alere Determine™ HIV-1/2 Ag/Ab Combo	Sensitivity
HIV-1	885	885	100.00 %
HIV-1 group O	7	7	100.00 %
HIV-1 non B subtypes*	154	154	100.00 %
HIV-2	133	133	100.00 %
Total	1179	1179	100.00 %

* Subtypes: A, C, D, F, G, H, J, K, and CRF AE, AG, AD, BD, 06, 09 and 011.

A total of 1179 confirmed HIV Ab positive specimens were tested (Table I). The diagnostic sensitivity of **Alere Determine™ HIV-1/2 Ag/Ab Combo** on this population of specimens is calculated to be 100%.

Table II
Comparison of results obtained with Alere Determine™ HIV-1/2 Ag/Ab Combo using matched specimens of whole blood (venipuncture and fingerstick), serum and plasma.

	Specimens				Correlation between matrices
	Type of Specimen				
No. of matched specimens tested	Serum	Plasma	Whole Blood (venipuncture)	Whole Blood (fingerstick)	
91	91	91	91	-	100 %
20	-	-	20	20	100 %
9	-	9	9	-	100 %
22	22	22	-	22	100 %
		Total Number of specimens tested			
142	113	122	120	42	100 %

Multiple (matched) specimens:

Seropositive specimens from a total of 142 individuals from Africa, Europe and South America were tested. Multiple (matched) specimens were obtained from several of these donors. From these 142 individuals, 113 serum specimens, 122 plasma specimens, 120 whole blood (venipuncture) specimens and 42 whole blood (fingerstick) specimens were obtained in various combinations.

Whole Blood (venipuncture) specimens:

120 whole blood (venipuncture) specimens were tested. Of these, 91 were matched with serum and plasma, 9 were matched pairs with plasma and 20 were matched pairs with whole blood (fingerstick) specimens.

Whole Blood (fingerstick) specimens:

42 whole blood (fingerstick) specimens were tested. Of these, 22 were matched with serum and plasma and 20 were matched pairs with whole blood (venipuncture) specimens.

The results obtained from all specimen matrices showed 100% correlation, demonstrating that **Alere Determine™ HIV-1/2 Ag/Ab Combo** gives identical results for these types of specimen matrices.

2. HIV seroconversion specimens:**Table III**
Seroconversion panels

Alere Determine™ HIV-1/2 Ag/Ab Combo results compared to Alere Determine™ HIV-1/2 Ag/Ab Combo			
33 seroconversion panels	Earlier detection (at least one bleed)	Equivalent detection (Same sample recognized as positive)	Later detection
Number of seroconversion panels			
23		10	0

A total of 33 seroconversion panels were studied and the **Alere Determine™ HIV-1/2 Ag/Ab Combo** results were compared to the results of the CE marked **Alere Determine™ HIV-1/2** (Table III).

With the exception of one panel, there was at least one specimen and up to 5 specimens reactive for the Ag bar, and recognized as acute infection (verified as positive for Ag according to the panel data sheet). **Alere Determine™ HIV-1/2 Ag/Ab Combo** detected HIV infection 2-20 days earlier than the **Alere Determine™ HIV-1/2** (3rd generation) antibody test, depending on the panel tested.

3. Early seroconversion:**Table IV**
Early seroconversion samples

Number of 42 early seroconversion samples	Alere Determine™ HIV-1/2 Ag/Ab Combo	Alere Determine™ HIV-1/2 (Ab detection)
	Positive	Positive
42	41	3

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


1. Pilcher C, Eron JJ, Galvin S, Gay S and Cohen MS (2004) Acute HIV revisited: new opportunities for treatment and prevention. *The Journal of Clinical Investigations.* 113(7): 937-945.
2. Respass RA, Rayfield MA and Dondero TJ (2001) Laboratory testing and rapid HIV assays: applications for HIV surveillance in hard-to-reach populations. *AIDS.* 15 Supplement 3: S49-S59.
3. Louie B, Wong E, Klausner JD, Liska S, Hecht F, Dowling T, Obeso M, Phillips SS and Pandori MW (2008) Assessment of rapid tests for detection of Human Immunodeficiency Virus-specific antibodies in recently infected individuals. *Journal of Clinical Microbiology.* 46(4): 1494-1497.
4. Lyons MS, Lindsell CJ, Hawkins DA, Raab DL, Trott AT and Fichtenbaum CJ (2008) Contributions to early HIV diagnosis among patients linked to care vary by testing venue. *BMC Public Health* 8: 220
5. Gürler LG, Hauser PH, Eberle J, et al. A New Subtype of Human Immunodeficiency Virus Type 1 (MVP-5180) from Cameroon. *Journal of Virology.* 1994; 68(3): 1581-1585.
6. Clinical and Laboratory Standards Institute. Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard-5th H4-A5 Vol. 24 No. 21
7. National Committee for Clinical Standards Clinical Waste Management: Approved Guideline. NCCLS Document GP5-A. Villanova, PA: NCCLS, 1993; 13(22): 1-18, 29-42.
8. US Environmental Protection Agency EPA Guide for Infectious Waste Management: Publication No. EPA/530-SW-86-014. Washington, DC: US Environmental Protection Agency, 1986: 1-1-5, R1-R3, A1-A24.
9. CDC: 2008 Case Definition; Human Immunodeficiency Virus Infection.
10. Delaney KP, Branson BM, et al. Evaluation of the Performance Characteristics of 6 Rapid HIV Antibody Tests. *Clinical Infectious Diseases.* 2011; 52(2): 257-263.
11. O'Connell RJ, Merritt TM, Malia JA, et al. Performance of the OraQuick Rapid Antibody Test for Diagnosis of Human Immunodeficiency Virus Type 1 Infection in Patients with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology.* 2003; 41(5):2153-2155.
12. O'Connell RJ, Agan BK, Anderson SA, Malia JA, Michael NL. Sensitivity of the Multispot HIV-1/HIV-2 Rapid Test Using Samples from Human Immunodeficiency Virus Type 1-Positive Individuals with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology.* 2006; 44(5): 1831-1833.

Advice Line

For further information, please contact your distributor, or call to one of the following Alere Product Support Care Centers:

Region	Phone	E-Mail Address
Europe & Middle East	+44 (0) 161 483 9032	EMEproductsupport@alere.com
Asia Pacific	+ (61) 7 3363 7711	APproductsupport@alere.com
Africa, Russia & CIS	+ (972) 8 9429 683	ARCIProductsupport@alere.com
Latin America	+ (57) 2 661 8797	Laproductsupport@alere.com

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