
BinaxNOW[®]







*Filaria*sis

Product Instructions
Consignes sur le produit
Produktanleitung
Istruzioni per l'uso del prodotto
Instruções sobre o produto
Instrucciones de uso del producto
Produktvejledning
Gebruiksaanwijzing voor het product
Produktanvisningar
Οδηγίες προϊόντος
Rimelige kompliserte produktinstrukser
تعليمات منتج
产品指示
제품 사용지침서

MATERIALS PROVIDED / MATÉRIAUX SUPPLÉER /
 GELIEFERTE MATERIALIEN / MATERIALE FORNITO /
 MATERIAIS FORNECIDOS / MATERIALES SUMINISTRADOS /
 MEDFØLGENE MATERIALER /
 MATERIEEL VOORZIEN / MATERIALEN FÖRSYNT /
 ARBEIDSMATERIALE FORSYNT /
 ΥΛΙΚΑ ΠΟΥ ΠΑΡΕΧΟΝΤΑΙ /
 محتويات علبة مجموعة الفحص
 测试包组件/유의 및 경고 사항

TEST PROCEDURE / ÉPREUVE PROCÉDURE / TESTVERFAHREN / ESAME PROCEDURA /
 PROCEDIMIENTO DE TESTE / PROCEDIMIENTO DE ENSAYO / TESTPROCEDURE / TOETS WIJZE VAN HANDELEN /
 PROV PROCEDUR / TEST FREMGANGSMÅTE / ΔΙΑΔΙΚΑΣΙΑ ΔΟΚΙΜΗΣ
 طريقة إجراء الفحص / 测试步骤 / 테스트 절차

Read result in window.
 Lire le résultat dans la fenêtre de lecture.
 Ergebnis im Fenster ablesen.
 Leggere il risultato nella finestra.
 Ver o resultado na janela.
 Vea el resultado en la ventana.
 Afleas resultatet i vinduet.
 Lees het resultaat in het venster af.
 Avläsning av resultat i fönstret.
 Les av resultatet i vinduet.
 Διαβάστε το αποτέλεσμα στο παράθυρο.

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ENGLISH

A rapid Immunochromatographic test (ICT) for the qualitative detection of *Wuchereria bancrofti* antigen.

For *in vitro* diagnostic use.

TEST PRINCIPLE

The BinaxNOW® *Filariasis* is an *in vitro* immunodiagnostic test for the detection of *Wuchereria bancrofti* antigen in whole blood, serum or plasma. The test utilizes an antibody specific for *W. bancrofti*. The antibody is attached to colloidal gold and impregnated into the pink and white sample pad. In the first step the sample is added to the top of the pink & white pad. If whole blood is used, the majority of the red blood cells will be retained in the white area of the pad and serum will flow forward into the pink area allowing any *W. bancrofti* antigen present to bind to the colloidal gold-labeled antibody. The antibody is immobilized in a line across the membrane. When the card is closed, the sample and the labeled antibody on the pink area of the pad contact the end of the membrane. The sample and the labeled antibody then migrate along the membrane crossing the immobilized antibody line. In a positive sample, any *W. bancrofti* antigen complexed with the gold-labeled antibody is captured by the antibody on the membrane and a pink line forms. In a negative sample, no gold-labeled antibody is captured by the antibody on the membrane and no pink line forms. A procedural control line will always appear in area C of the test window if the test has been performed correctly.

SPECIMEN COLLECTION

Capillary or venous blood collected in a heparin coated container may be used. Plasma (heparin) or serum may also be used in the test.

- Finger prick blood collection: Collect the sample in a 100µl calibrated capillary. When sampling, care should be taken to avoid air bubbles in the capillary used to collect the blood.
- Collect venous blood by the standard venipuncture procedure into a heparin tube.

SAMPLE PRESERVATION

Blood samples collected in a capillary should be used immediately. Venous blood can be used for up to three days if stored at 2° to 8°C. Serum or plasma can be stored for up to two weeks at 2° to 8°C.

PRECAUTIONS AND WARNINGS

- Optimal results will be obtained by strict adherence to this protocol. Do not reuse test cards.
- Biological contamination of dispensing equipment or containers can lead to false results. Observe established precautions against microbiological and serological hazards in specimen handling, disposal and throughout all procedures.
- Do not use components past expiration date. Keep storage boxes dry.

LIMITATIONS OF THE PROCEDURE

This test is structured to indicate the presence or absence of *W. bancrofti* antigen in the sample. The absence of antigen does not exclude *Filariasis* caused by other nematode species.

KIT CONTENTS

- 25 individually packaged test cards.
- 25 capillaries.

KIT STORAGE

Store kit at 2-30°C. The BinaxNOW® *Filariasis* test kit is stable until the expiration date marked on its outer packaging when stored as specified.

TEST PROCEDURE

- Remove the test card from the pouch just prior to use. Open the card and lay it flat on the work surface.
- To ensure good blood flow and performance, **completely** fill the capillary tube with blood.
Note: When using venous blood, serum, or plasma, either a calibrated capillary or calibrated pipette may be used.
- Add the 100µl of sample **slowly** from the capillary or pipette onto the top of the pink & white pad as shown. **Note: See arrow on test device showing correct placement of sample at TOP of pad.** Should there be sample remaining in the capillary that will not flow out freely, the tip may be gently pressed against the pad. **IMPORTANT: Wait for each drop to soak in before adding the next drop. Incorrect addition of sample may result in device failure.**
- Wait until the sample has flowed into the pink area and it is completely wet (this should take about 30 seconds to 1 minute).
- Remove and discard the adhesive liner. Ensure that the adhesive on the right hand side of the test card is exposed.
- Close the card. To ensure good test flow, press very firmly along the entire area to the right of the window. Start timing.
- Read the result through the viewing window at 10 minutes. Refer to the details given in this insert for test interpretation.

RESULTS INTERPRETATION

Positive Test Result

The test is positive if two lines (T and C) are seen in the viewing window. Any pink line in the T area indicates a positive test result. The test is positive even when the T line appears lighter or darker than the C line.

Negative Test Result

The test is negative if only the C line is seen. To ensure that low positive samples have had sufficient time to develop, a negative result should not be recorded until 10 minutes have elapsed from when the card is closed.

Invalid Test Result

The test is invalid if the C line does not appear. If this occurs, the test should be repeated.

ORDERING and CONTACT INFORMATION

Reorder numbers:

#620-000: BinaxNOW® Filariasis Test Kit

OUS +1-321-441-7200

Technical Support

Advice Line

Further information can be obtained from your distributor, or by contacting Alere™ Technical Support on:

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© 2014 Alere. All rights reserved.
The Alere Logo, Alere and BinaxNOW are trademarks of the Alere group of companies.
All trademarks referenced are trademarks of the respective owners.
Alere Scarborough, Inc.