


# E. HISTOLYTICA QUIK CHEK™ Test Procedure


## 1 Add to a test tube:

- 500  $\mu$ L *Diluent*
- 1 drop *Conjugate*
- 25  $\mu$ L (1<sup>st</sup> graduation) of specimen\*
- Mix thoroughly




\*For specimens in transport media, use 400  $\mu$ L (1<sup>st</sup> graduation) of *Diluent* and 100  $\mu$ L (2<sup>nd</sup> graduation on pipette) specimen.


## 2 Transfer 500 $\mu$ L (highest graduation) from tube to small *Sample Well*. Keep the cassette at room temperature and wait 15 minutes.



## 3 Add 300 $\mu$ L *Wash Buffer* to large *Reaction Window*. Allow to completely absorb.



## 4 Add 2 drops *Substrate* to large *Reaction Window*. Keep the cassette at room temperature and wait 10 minutes. Read results.



## Interpretation

A positive result may be interpreted at any time during the 10 minute incubation.

A test cannot be ruled negative or invalid until 10 minutes have passed.



**POSITIVE**



**NEGATIVE**



**INVALID**



**INVALID**

If there is no line at C (control) the test is invalid. See package insert for additional interpretation of results.

Distributed by:



Developed and Manufactured by:



This chart does not contain the complete instructions for use of the *E. HISTOLYTICA QUIK CHEK* test. For proper use of the assay, please read the package insert.

Made in the USA. The Alere Logo and Alere are trademarks of the Alere group of companies. *E. HISTOLYTICA QUIK CHEK*, the TECHLAB Logo and TECHLAB are trademarks of TECHLAB, Inc., under license. © 2017 TECHLAB, Inc. All rights reserved.

Part No. 91-409-02  
Issued: 7/2017