AFINION™ CRP

EN 5-13

EL 23-32

FR 43-52

NL 62-70

DE 14-22

ES 33-42

IT 53-61

PT 71-80

REF 1116787

1116908 Rev. B 2020/08

SYMBOLS/ SYMBOLE/ ΣΥΜΒΟΛΑ/ SÍMBOLOS/ SYMBOLES/ SIMBOLI/ SYMBOLEN/ SÍMBOLOS

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
EN AFINION™ CRP

For use with the Alere Afinion™ AS100 Analyzer/Afinion™ 2 Analyzer. For professional near-patient testing and laboratory use.

PRODUCT DESCRIPTION

Intended use

Afinion™ CRP is an *in vitro* diagnostic test for the quantitative determination of C-reactive protein (CRP) in human whole blood and in human serum and plasma. The measurement of CRP provides information for the detection and evaluation of infection, tissue injury, inflammatory disorders and associated diseases.

Summary and explanation of the test

C-reactive protein (CRP) is one of the cytokine induced acute-phase proteins, the levels of which rise during a general, unspecific response to infections and non-infectious inflammatory processes^{5,7}. In healthy persons the serum or plasma CRP levels are below 5 mg/L¹⁻³. This threshold is often exceeded within four to eight hours after an acute inflammatory event, with CRP values reaching approximately 20 to 500 mg/L⁴. As elevated CRP levels are always associated with pathological changes, the CRP assay provides information for the diagnosis, therapy, and monitoring of inflammatory diseases⁴⁻⁷.

Principle of the assay

Afinion CRP is a solid phase immunochemical assay. The Afinion CRP Test Cartridge contains all the reagents necessary for measuring the CRP concentration in whole blood, serum or plasma. The sample material is collected using the integrated sampling device and the test cartridge is placed in the Afinion Analyzer. The sample is then automatically diluted with a liquid that also lyses the blood cells. The sample mixture is passed through a membrane coated with anti-CRP antibodies; the CRP in the sample is concentrated onto this membrane. A solution containing anti-CRP antibodies conjugated with ultra-small gold particles is passed through the membrane. The gold-antibody conjugate binds to the immobilized CRP on the membrane, which turns red-brown. Excess gold-antibody conjugate is removed by a washing solution. The Afinion Analyzer measures the colour intensity of the membrane, and this is proportional to the amount of CRP in the sample. The CRP concentration is displayed on the analyzer screen.

Kit contents (per 15 tests unit)


- 15 Test cartridges packaged separately in foil pouches with a desiccant bag
- 1 Package insert

Materials required but not provided with the kit

- Alere Afinion AS100 Analyzer/ Afinion 2 Analyzer
- Afinion CRP Control
- Standard blood collection equipment

Description of the test cartridge

The main components of the test cartridge are the sampling device and the reagent container. The test cartridge has a handle, a barcode label with lot-specific information and an area for sample ID. See figure and table below.

Component	Function/composition
1 Sampling device a. Closed position b. Lifted position	For collection of patient sample or control.
2 Capillary	2.5 µL glass capillary to be filled with sample material.
3 Reaction wells a. Capillary wiper b. Membrane tube c. Conjugate solution d. Dilution liquid e. Washing solution	Contain reagents necessary for one test: Plastic laminated paper. Tube with a nitrocellulose membrane coated with anti-CRP antibodies. Anti-CRP antibodies labelled with ultra-small gold particles. TAPS buffer solution. TAPS buffer solution and detergents.
4 Handle	For correct finger grip.
5 Barcode label	Contains assay- and lot-specific information for the analyzer.
6 Optical reading area	Area for transmission measurement.
7 ID area	Space for written or labelled sample identification.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use.
- Do not use test cartridges after the expiry date or if the test cartridges have not been stored in accordance with recommendations.
- Do not use the test cartridge if the foil pouch or the test cartridge has been damaged.
- Each foil pouch contains a desiccant bag with 1 g silica gel. This material shall not be used in the assay. Discard the desiccant bag in a suitable container. Do not swallow.
- Do not use the test cartridge if the desiccant bag is damaged and desiccant particles are found on the test cartridge.
- In case of leakage from the test cartridge, avoid contact with eyes and skin. Wash with plenty of water.
- Do not reuse any part of the test cartridge.

- The used test cartridges, sampling equipment, patient samples and controls are potentially infectious. The test cartridges should be disposed of immediately after use. Proper handling and disposal methods should be followed in accordance with local or national regulations. Use gloves.

STORAGE

Refrigerated storage (2-8°C)

- The Afinion CRP Test Cartridges are stable until the expiry date only when stored refrigerated in sealed foil pouches.
- Do not freeze.

Room temperature storage (15-25°C)

- The Afinion CRP Test Cartridges can be stored in unopened foil pouches at room temperature for 4 weeks. Note the date of removal from the refrigerator on the kit container.
- Avoid exposure to direct sunlight.

Opened foil pouch

- The test cartridge must be used within 10 minutes after opening the foil pouch.
- Avoid exposure to direct sunlight.
- Avoid relative humidity above 90 %.

SAMPLE MATERIAL

The following sample materials can be used with the Afinion CRP Test:

- Capillary blood
- Venous whole blood with anticoagulants (K₂-EDTA, Li-Heparin, Na-Heparin)
- Serum
- Plasma (K₂-EDTA, Li-Heparin, Na-Heparin)
- Afinion CRP Control

Sample storage

- Capillary blood without anticoagulants cannot be stored.
- Venous whole blood with anticoagulants (K₂-EDTA, Li-Heparin, Na-Heparin) can be stored refrigerated for 3 days. Do not freeze.
- Serum and plasma can be stored refrigerated for 10 days.
- Serum and plasma can be stored frozen for one year if tubes are sealed properly.
- Consult the Afinion CRP Control Package Insert for storage of the control materials.

PREPARING FOR ANALYSIS

Consult the Afinion Analyzer User Manual for instructions on how to analyse a patient or control sample. The Afinion CRP Quick Guide also provides an illustrated step-by-step procedure.

- The Afinion CRP Test Cartridge must reach an operating temperature of 15-30°C before use. Upon removal from refrigerated storage, leave the test cartridge in unopened foil pouch for at least 15 minutes.
- Open the foil pouch just before use.
- Do not touch the test cartridge optical reading area.
- Label the test cartridge with patient or control ID. Use the dedicated ID area.

COLLECTING A SAMPLE

Capillary blood

- Make sure that the finger is clean, warm and dry.
- Use a suitable lancet to prick the finger.
- Allow a good drop of blood to form before sampling.

Venous blood, serum or plasma




- Patient samples stored refrigerated can be used without equilibration to room temperature.
- Mix the sample material well by inverting the tube 8-10 times before collecting a sample.

AFINION™ CRP Control

- The control material can be used without equilibration to room temperature.
- Mix the control material well by inverting the vial 8-10 times before collecting a sample.

Filling the capillary

- Remove the sampling device from the test cartridge.
- Fill the capillary; bring the tip of the capillary just beneath the surface of the patient sample (a, b) or control material (c). Be sure that the capillary is completely filled, see arrow (d). It is not possible to overfill. Avoid air bubbles and excess sample on the outside of the capillary. Do not wipe off the capillary.
- Reinsert the sampling device immediately into the test cartridge.
- Once the capillary is filled with sample, the analysis of the test cartridge must start within 1 minute. If the test cartridge is stored too long before analysis, the sample material may dry or coagulate. An information code will then be displayed. Consult the Afinion Analyzer User Manual.
- Do not use a test cartridge that has been accidentally dropped on the floor or lab bench after sample collection.

ANALYSING A SAMPLE

- Analyse the test cartridge following the procedure described in the Afinion Analyzer User Manual.
- The analysis time is 3-4 minutes.

INTERPRETATION OF RESULTS

Interpret the Afinion CRP test results with careful consideration to the patient's medical history, clinical examinations and other laboratory results. If the test result is questionable or if clinical signs and symptoms appear inconsistent with the test result, analyze the Afinion CRP Controls and retest the sample using a new Afinion CRP Test Cartridge. If the result is still questionable, confirm the result with another method. Analyse control materials frequently to verify the performance of the Afinion Analyzer System.

Measuring range

- Whole blood samples:
CRP 5-200 mg/L.
CRP < 5 mg/L is displayed if the CRP concentration is below 5 mg/L.
CRP > 200 mg/L is displayed if the CRP concentration is above 200 mg/L.
- Serum and plasma samples:
CRP 5-160 mg/L.
CRP < 5 mg/L is displayed if the CRP concentration is below 5 mg/L.
CRP > 160 mg/L is displayed if the CRP concentration is above 160 mg/L.

If the CRP concentration is above the measuring range, serum or plasma samples can be diluted with 0.9 % saline up to 4 times (1:4) and retested. Diluted whole blood cannot be used.

Standardisation

Afinion CRP is calibrated against the ERM®-DA/IFCC reference preparation.

Reference range

CRP < 5 mg/L¹⁻³

Automatic correction for hematocrit (Hct)

When whole blood is used as sample material, the measured CRP value is automatically corrected according to the sample's Hct level (limited to range 20-60 %), and the result is displayed as the plasma-CRP concentration. A CRP result obtained from analysis of a whole blood sample is thus directly comparable to the CRP concentration measured in the serum or plasma fraction of the same sample. If the Hct value is outside the range 20-60 %, no CRP test result will be reported and an information code will be displayed (see "Troubleshooting"). In these cases serum or plasma samples are recommended for CRP analysis.

Analytical specificity

Monoclonal antibodies specific to human CRP are used in Afinion CRP. Serum Amyloid P component (SAP), a protein with close homology to CRP, does not cross-react with CRP in the Afinion CRP Test System.

Interference

- No significant interference (0 ± 4 mg/L at CRP < 40 mg/L, 0 ± 10 % at CRP ≥ 40 mg/L) was observed up to the following concentrations:

Bilirubin	502 $\mu\text{mol/L}$ (29.3 mg/dL)
Cholesterol	10 mmol/L (400 mg/dL)
Leukocytes*	$50 \times 10^9/\text{L}$
Rheumatoid factor (RF)	1003 IU/mL
Serum Amyloid P component	100 mg/L
Triglycerides	8 mmol/L (700 mg/dL)

Over-the-counter and prescription drugs:

Acetaminophen (Paracetamol)	1334 $\mu\text{mol/L}$ (0.2 mg/mL)
Acetylsalicylic acid	3611 $\mu\text{mol/L}$ (0.7 mg/mL)
Amoxicillin (Penicillins)	208 $\mu\text{mol/L}$ (0.09 mg/mL)
Cefotaxime (Cephalosporines)	673 $\mu\text{mol/L}$ (0.32 mg/mL)
Ciprofloxacin (Quinolones)	32 $\mu\text{mol/L}$ (0.01 mg/mL)
Clindamycin (Lincomycines)	90 $\mu\text{mol/L}$ (0.05 mg/mL)
Erythromycin (Macrolides)	83 $\mu\text{mol/L}$ (0.06 mg/mL)
Ibuprofen	2425 $\mu\text{mol/L}$ (0.5 mg/mL)
Salicylic acid	4347 $\mu\text{mol/L}$ (0.6 mg/mL)
Sulfamethoxazole (Sulfonamides)	1589 $\mu\text{mol/L}$ (0.4 mg/mL)
Tetracyclin (Tetracyclines)	34 $\mu\text{mol/L}$ (0.02 mg/mL)

* For neonates no interference from leukocytes was observed up to $30 \times 10^9/\text{L}$. A positive interference has been observed above this level of leukocytes.

- No “Hook” effect is observed at CRP concentrations up to 2000 mg/L.
- Anticoagulants (K_2 -EDTA, Li-Heparin, Na-Heparin) at concentrations normally used in blood collection tubes do not interfere.

Important! It is possible that other substances and/or factors not listed above may interfere with the test and cause false results.

Limitations of the test

- Diluted whole blood samples cannot be used with Afinion CRP. If the CRP concentration is above the measuring range, diluted serum or plasma can be used.
- If the Hct value is outside the range 20-60 %, no CRP test result will be reported and an information code will be displayed (see “Troubleshooting”). In these cases serum or plasma samples are recommended for CRP analysis.
- Afinion CRP is not intended for measuring CRP as a risk marker of coronary heart disease. For this purpose, a *high sensitivity* CRP test has to be used.

QUALITY CONTROL

Quality control testing should be done to confirm that your Afinion Analyzer System is working properly and providing reliable results. Only when controls are used routinely and the values are within acceptable ranges can accurate results be assured for patient samples.

It is recommended to keep a permanent record of all quality control results. The Afinion Analyzer automatically stores the control results in a separate log. Consult the Afinion Analyzer User Manual.

Choosing control material

Afinion CRP Control from Abbott is recommended for routine quality control testing. Consult the Afinion CRP Control Package Insert.

If controls from another supplier are used, the precision must be determined and acceptable ranges for the Afinion Analyzer System must be established.

Frequency of control testing

Controls should be analysed:

- anytime an unexpected test result is obtained.
- with each shipment of Afinion CRP Test Kits.
- with each new lot of Afinion CRP Test Kits.
- when training new operators in correct use of the Afinion CRP and the Afinion Analyzer.
- in compliance with national or local regulations.

Verifying the control results

The measured value should be within the acceptable limits stated for the control material. Consult the Afinion CRP Control Package Insert.

If the result obtained for the Afinion CRP Control is outside the acceptable limits, make sure that:

- patient samples are not analysed until control results are within acceptable limits.
- the control vial has not passed its expiry date.
- the control vial has not been used for more than 4 weeks.
- the control vial and Afinion CRP Test Cartridge have been stored according to recommendations.
- there is no evidence of bacterial or fungal contamination of the control vial.

Correct any procedural error and retest the control material.

If no procedural errors are detected:

- Examine the laboratory's quality control record to investigate the frequency of control failures.
- Ensure that there is no trend in out-of-range quality control results.
- Retest the control material using a new control vial.
- Patient results must be declared invalid when controls do not perform as expected. Contact your local supplier for advice before analysing patient samples.

TROUBLESHOOTING

To ensure that correct CRP results are reported, the Afinion Analyzer performs optical, electronic and mechanical controls of the capillary, the test cartridge and all individual processing steps during the course of each analysis. When problems are detected by the built-in failsafe mechanisms, the analyzer terminates the test and displays an information code.

The table below contains the Afinion CRP specific information codes. Consult the Afinion Analyzer User Manual for information codes not listed in this table.

Code #	Cause
101	Hematocrit below 20 %
102	Hematocrit above 60 %

PERFORMANCE CHARACTERISTICS

Method comparison

Afinion CRP used with Alere Afinion AS100 Analyzer was compared with an automated immunoturbidimetric laboratory method. Li-heparin plasma and EDTA whole blood analysed with Afinion CRP were compared with Li-heparin plasma analysed with the laboratory method. The results are presented in Table 1.

Table 1: Comparison between Afinion CRP and a laboratory method (Lab) using Li-heparin plasma (pl) and EDTA whole blood (wb).

AFINION™ vs. laboratory method	y: AFINION™ (pl) x: Lab (pl)	y: AFINION™ (wb) x: Lab (pl)
Number of samples	102	110
Regression line	$y = 0.95x - 0.38$	$y = 1.00x - 0.02$
Correlation coefficient	$r = 1.00$	$r = 1.00$

Precision

Precision studies were conducted according to CLSI Guideline EP05-A3. Within-device precision was determined for 4 plasma samples, 4 serum samples and 2 control samples assayed for 20 days, each sample was assayed in duplicate, twice a day. Four whole blood samples were assayed for 3 days, each sample was assayed in 28 replicates per day. The data are summarized in Table 2.

Table 2: Within-device precision. SD=Standard deviation, CV=Coefficient of Variation.

Sample	Number of days	Mean CRP (mg/L)	Within-device precision	
			SD (mg/L)	CV (%)
Whole blood 1	3	9.1	0.4	4.5
Whole blood 2	3	19.2	0.8	4.2
Whole blood 3	3	48.5	1.7	3.5
Whole blood 4	3	174.8	6.3	3.6
Serum 1	20	7.4	0.3	3.5
Serum 2	20	18.9	0.6	3.1
Serum 3	20	47.7	1.1	2.4
Serum 4	20	148.1	5.2	3.5
Plasma 1	20	7.5	0.2	3.0
Plasma 2	20	17.8	0.6	3.1
Plasma 3	20	48.4	1.4	2.8
Plasma 4	20	141.5	4.4	3.1
Afinion CRP Control C I	20	20.2	0.6	2.9
Afinion CRP Control C II	20	57.6	1.4	2.5

Performance testing with the AFINION™ 2 Analyzer

The performance of Afinion CRP obtained with the Afinion 2 Analyzer has been demonstrated to be equivalent to the performance obtained with the Alere Afinion AS100 Analyzer.

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1. Claus DR, Osmand AP, Gewurz H. Radioimmunoassay of human C-reactive protein and levels in normal sera. *J Lab Clin Med* 1976; 87:120-8.
2. Dati F, Schumann G et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference ranges for 14 proteins in serum based on the standardization against the IFCC/BCR/CAP Reference Material (CRM 470). International Federation of Clinical Chemistry. Community Bureau of Reference of the Commission of the European Communities. College of American Pathologists.
Eur J Clin Chem Clin Biochem 1996; 34(6):517-20.
3. Dati F, Johnson AM, Whicher JT. The Existing Interim Consensus Reference Ranges and the Future Approach. *Clin Chem Lab Med* 2001; 39(11):1134-6.
4. Morley J, Kushner I. Serum C-reactive protein levels in disease. *Ann NY Acad Sci* 1982; 389:406-18.
5. Kushner I, Rzewnicki DL. The acute phase response: general aspects. *Baillieres Clin Rheumatol.* 1994 Aug;8(3):513-30.
6. Pepys MB, Hirschfield GM. C-reactive protein: a critical update. *J Clin Invest.* 2003 Jun;111(12):1805-12.
7. van Leeuwen MA, van Rijswijk MH. Acute phase proteins in the monitoring of inflammatory disorders. *Baillieres Clin Rheumatol.* 1994 Aug;8(3):531-52.

 **Abbott Diagnostics Technologies AS**
Kjelsåsveien 161
P.O. Box 6863 Rodeløkka
NO-0504 Oslo, Norway
www.abbott.com/poct

ISO 13485 certified company

© 2020 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.

1116908 Rev. B 2020/08