AFINION™ HbA1c

EN 5-14

EL 25-34

FR 45-54

NL 65-74

DE 15-24

ES 35-44

IT 55-64

PT 75-84

REF 1116795, 1117196

1117194 Rev. A 2021/10

SYMBOLS /SYMBOLE/ ΣΥΜΒΟΛΑ/ SÍMBOLOS/ SYMBOLES/SIMBOLI/ SYMBOLEN/SÍMBOLOS

CE

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EN AFINION™ HbA1c

For use with the Alere Afinion™ AS100 Analyzer/Afinion™ 2 Analyzer.
For professional near-patient testing and laboratory use.

PRODUCT DESCRIPTION

Intended use

Afinion™ HbA1c is an *in vitro* diagnostic test for quantitative determination of glycated hemoglobin (hemoglobin A1c, HbA1c) in human whole blood. The measure of HbA1c is recommended as a marker of long-term metabolic control in persons with diabetes mellitus. This test can be used as an aid in the diagnosis of diabetes and as an aid in identifying patients who may be at risk for developing diabetes.

Summary and explanation of the test

The human erythrocyte is freely permeable to glucose. Within each erythrocyte a slow, continuous, non-enzymatic process between hemoglobin A and various sugars takes place. The product formed is known as glycated hemoglobin, or glycohemoglobin¹. The chronic elevated blood sugar level of persons with diabetes mellitus will over time cause damage to the small vessels of the body. This damage develops slowly over years and is known to cause late complications. Good metabolic control, i.e. lowering the HbA1c concentration, has proven to delay the onset and slow the progression of diabetes late complications^{3,4,5}.

An International Expert Committee has concluded that measurements of HbA1c can be used to diagnose diabetes mellitus. When in agreement with national regulations, Afinion HbA1c can be used as an aid in the diagnosis of diabetes and as an aid in identifying patients who may be at risk for developing diabetes².

Principle of the assay

Afinion HbA1c is a fully automated boronate affinity assay for the determination of hemoglobin A1c in human whole blood.

The Afinion HbA1c Test Cartridge contains all the reagents necessary for the determination of the HbA1c concentration. The sample material is collected using the integrated sampling device and the test cartridge is placed in the Afinion Analyzer. The blood sample is then automatically diluted and mixed with a liquid that releases hemoglobin from the erythrocytes. The hemoglobin precipitates. This sample mixture is transferred to a blue boronic acid conjugate, which binds to the cis-diols of glycated hemoglobin. This reaction mixture is soaked through a filter membrane and all precipitated hemoglobin, conjugate-bound and unbound (i.e. glycated and non-glycated hemoglobin) remains on the membrane. Any excess of conjugate is removed with a washing reagent.

The Afinion Analyzer evaluates the precipitate on the membrane. By measuring the reflectance, the blue (glycated hemoglobin) and the red (total hemoglobin) colour intensities are evaluated, the ratio between them being proportional to the percentage of HbA1c in the sample. The HbA1c concentration is displayed on the Afinion Analyzer in the units mmol/mol, percentage (%), as the calculated estimated average glucose (eAG) or in combinations of these.

Kit contents (per 15 tests unit)


- 15 Test cartridges packaged separately in foil pouches with a desiccant bag
- 1 Package insert

Materials required but not provided with the kit

- Alere Afinion AS100 Analyzer/Afinion 2 Analyzer
- Afinion HbA1c Control
- Standard blood collection equipment

Description of the test cartridge

The main components of the test cartridge are the sampling device and the reagent container. The test cartridge has a handle, a barcode label with lot-specific information and an area for sample ID. See figure and table below.

Component	Function/composition
1 Sampling device a. Closed position b. Lifted position	For collection of patient sample or control.
2 Capillary	1.5 µL capillary to be filled with sample material.
3 Reaction wells a. Conjugate b. Membrane tube c. Washing solution d. Reconstitution reagent e. Empty	Contains reagents necessary for one test: Blue boronic acid conjugate. Tube with a polyethersulfone membrane. Morpholine buffered sodium chloride with detergents and preservative. HEPES buffered sodium chloride with lysis and precipitation agents. N/A
4 Handle	For correct finger grip.
5 Barcode label	Contains assay- and lot-specific information for the analyzer.
6 Optical reading area	Area for transmission measurement.
7 ID area	Space for written or labelled sample identification.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use.
- Do not use test cartridges after the expiry date or if the test cartridges have not been stored in accordance with recommendations.
- Do not use the test cartridge if the foil pouch or the test cartridge has been damaged.
- Each foil pouch contains a desiccant bag with 1 g silica gel. This material shall not be used in the assay. Discard the desiccant bag in a suitable container. Do not swallow.
- Do not use the test cartridge if the desiccant bag is damaged and desiccant particles are found on the test cartridge.

- Do not reuse any part of the test cartridge.
- The test cartridge contains sodium azide (< 0.1 %) as a preservative. In case of leakage, avoid contact with eyes and skin. Wash with plenty of water.
- The used test cartridges, sampling equipment, patient samples and controls are potentially infectious. The test cartridges should be disposed of immediately after use. Proper handling and disposal methods should be followed in accordance with local or national regulations. Please also refer to the Safety Data Sheet available at <https://www.globalpointofcare.abbott>.
- Use gloves.

STORAGE

Refrigerated storage (2-8°C)

- The Afinion HbA1c Test Cartridges are stable until the expiry date only when stored refrigerated in sealed foil pouches.
- Do not freeze.

Room temperature storage (15-25°C)

- The Afinion HbA1c Test Cartridges can be stored in unopened foil pouches at room temperature for 90 days. Note the date of removal from the refrigerator on the kit container.
- Avoid exposure to direct sunlight.

Opened foil pouch

- The test cartridge must be used within 10 minutes after opening the foil pouch.
- Avoid exposure to direct sunlight.
- Avoid relative humidity above 90 % (non-condensing).

SAMPLE MATERIAL

The following sample materials can be used with the Afinion HbA1c Test:

- Capillary blood (from finger prick)
- Venous whole blood with anticoagulants (K₂-EDTA, K₃-EDTA, Li-Heparin, Na-Heparin, NaF/Na₂-EDTA, NaF/K-oxalate, Na-citrate)
- Afinion HbA1c Control

Sample storage

- Capillary blood samples cannot be stored.
- Venous whole blood with anticoagulants (K₂-EDTA, K₃-EDTA, Li-heparin, Na-Heparin, Na-citrate) can be stored refrigerated (2-8°C) for 10 days or at room temperature (18-30°C) for 8 hours. Do not freeze.
- Consult the Afinion HbA1c Control Package Insert for storage of control materials.

PREPARING FOR ANALYSIS

Consult the Afinion Analyzer User Manual for detailed instructions on how to analyse a patient sample or control. The Afinion HbA1c Quick Guide also provides an illustrated step-by-step procedure.

- The Afinion HbA1c Test Cartridge must reach an operating temperature of 18-30°C before use. Upon removal from refrigerated storage, leave the test cartridge in unopened foil pouch for at least 15 minutes.
- Open the foil pouch just before use.
- Do not touch the test cartridge optical reading area.
- Label the test cartridge with patient or control ID. Use the dedicated ID area.

COLLECTING A SAMPLE

Capillary blood

- Make sure that the finger is clean, warm and dry.
- Use a suitable lancet to prick the finger.
- Allow a good drop of blood to form before sampling.

Venous blood

- Patient samples stored refrigerated can be used without equilibration to room temperature.
- Mix the sample material well by inverting the tube 8-10 times before collecting a sample.






AFINION™ HbA1c Control

- Allow the control material to reach room temperature (15-25°C) before use, which takes approximately 45 minutes.
- Mix the control material thoroughly by shaking the vial for 30 seconds. A whirl mixer may be used.
- The sample can be extracted from the vial or the cap.

Filling the capillary

The Afinion HbA1c Test Cartridge is equipped with either of the two models A or B of the sampling device (a). The two models have slightly different shape, but the filling volume of the capillary and the handling of the device is exactly the same:

- Remove the sampling device from the test cartridge.
- Fill the capillary; bring the tip of the capillary just beneath the surface of the patient sample (b, c) or control material (d). The capillary is automatically filled to the end, see arrow (e). It is not possible to overfill. Avoid air bubbles and excess sample on the outside of the capillary. Do not wipe off the capillary.
- Immediately replace the sampling device into the test cartridge.
- The analysis of the test cartridge should start within 2-3 minutes. If the test cartridge is stored too long before analysis, the sample material may dry or coagulate. An information code will then be displayed. Consult the Afinion Analyzer User Manual.
- Do not use a test cartridge that has been accidentally dropped on the floor or lab bench after sample collection.

ANALYSING A SAMPLE

- Analyse the test cartridge following the procedure described in the Afinion Analyzer User Manual.
- The analysis time is approximately 3.5 minutes.

INTERPRETATION OF RESULTS

Interpret the Afinion HbA1c test results with careful consideration to the patient’s medical history, clinical examinations and other laboratory results. If the test result is questionable or if clinical signs and symptoms appear inconsistent with the test result, analyze the Afinion HbA1c Controls and retest the sample using a new Afinion HbA1c Test Cartridge. If the result is still questionable, confirm the result with another method. Analyse control materials frequently to verify the performance of the Afinion System.

HbA1c measuring units

Three different measuring units are in use for reporting HbA1c test results⁸:

- **mmol/mol** – the HbA1c values are aligned to the IFCC Reference Method for Measurements of HbA1c^{6,7}
- **percentage (%)** – the HbA1c values are aligned to the assay used in the DCCT study, also known as NGSP-HbA1c^{3,8,9}
- **estimated average glucose (eAG)** – the HbA1c value is converted to an equivalent average glucose value¹⁰

Linear relationships have been established between the three units^{9,10}:

$$\begin{aligned} \text{NGSP-HbA1c (\%)} &= \text{DCCT-HbA1c (\%)} = 0.092 \times \text{IFCC-HbA1c (mmol/mol)} + 2.15 \\ \text{eAG (mmol/L)} &= 1.59 \times \text{DCCT-HbA1c (\%)} - 2.59 \end{aligned}$$

Measuring range

The Afinion Analyzer displays the HbA1c value/concentration in mmol/mol, percentage (%), as the calculated estimated average glucose (eAG) or in combinations of these:

	HbA1c	HbA1c	eAG
Unit	mmol/mol	%	mmol/L
Measuring range	20-140	4.0-15.0	3.8-21.3

Values outside the HbA1c measuring range

Valid for SW ≤ 7.03 (Afinion AS100) and SW ≤ 21.09 (Afinion 2):

If the patient’s HbA1c value is outside the measuring range, no test result will be reported and an information code will be displayed (see “Troubleshooting”).

Valid for SW ≥ 7.04 (Afinion AS100) and SW ≥ 21.10 (Afinion 2):

- Test result below the measuring range is displayed as less than (<) the lower limit.
- Test result above the measuring range is displayed as greater than (>) the upper limit.

 **Important!** Patient’s HbA1c values should be reported in units consistent with the national recommendations. Please contact your local supplier if the national recommendation is unknown. Consult the Afinion Analyzer User Manual for instructions on how to change the HbA1c measuring unit.

Standardisation

Afinion HbA1c is traceable to the IFCC Reference Method for Measurement HbA1c⁶.

Expected values

Patients with HbA1c levels between 5.7 and 6.4 % (39 and 47 mmol/mol) are identified as being at risk for developing diabetes. The diagnostic cut-off is 6.5 % (48 mmol/mol)².

Analytical specificity

Afinion HbA1c measures the total glycated hemoglobin and reports the HbA1c value. No significant interference (< 7%) was observed for samples with hemoglobin (Hb) variants and hemoglobin derivatives up to the following concentrations:

• HbA2	5.7 %
• HbAC	36 %
• HbAD	42 %
• HbAE	26 %
• HbAS	42 %
• HbF	10.4 %
• Acetylated Hb	4.6 mg/mL
• Carbamylated Hb	13.8 mg/mL
• Labile (pre-glycated) Hb	11.4 mg/mL

Interference

No significant interference (< 7%) was observed up to the following concentrations:

• Bilirubin conjugated	600 mg/L
• Bilirubin unconjugated	600 mg/L
• Glucose	10 g/L
• Lipids (as Intralipid)	10 g/L
• Rheumatoid factor	780 000 IU/L
• Total protein	150 g/L
• Glycated albumin	7.7 g/L

Over-the-counter and prescription drugs:

• Acetaminophen	200 mg/L
• Acetylcysteine	1663 mg/L
• Acetylsalicylic acid	1000 mg/L
• Ampicillin	1000 mg/L
• Ascorbic acid	300 mg/L
• Cefoxitin	2500 mg/L
• Cyclosporine A	5 mg/L
• Cyclosporine C	5 mg/L
• Doxycycline	50 mg/L
• Glyburide	1.9 mg/L
• Heparin	5000 U/L
• Ibuprofen	500 mg/L
• Levodopa	20 mg/L
• Metformin	40 mg/L
• Methyldopa	20 mg/L
• Metronidazole	200 mg/L
• Phenylbutazone	400 mg/L
• Rifampicin	64 mg/L
• Salicylic acid	599 mg/L
• Theophylline	100 mg/L
• Hemolysis (<i>in vitro</i>)	14 %
• Anticoagulants (K ₂ -EDTA, K ₃ -EDTA, Li-Heparin, Na-Heparin, NaF/Na ₂ -EDTA, NaF/K-oxalate and Na-citrate) at concentrations normally used in blood collection tubes do not interfere.	

Important! It is possible that other substances and/or factors not listed above may interfere with the test and cause false results.

Limitations of the test

- Any cause of shortened erythrocyte life span will reduce exposure of erythrocytes to glucose, resulting in a decrease in HbA1c values, regardless of the method used. Caution should be used when interpreting the HbA1c results from patients with conditions such as hemolytic anemia or other hemolytic diseases, homozygous sickle cell trait, pregnancy, blood loss, polycythemia, iron deficiency etc.
- This test should not be used to diagnose:
 - diabetes during pregnancy
 - patients with an elevated fetal hemoglobin (HbF >10 %) such as hereditary persistence of fetal hemoglobin (HPFH)
 - patients with a hemoglobinopathy but normal red cell turnover (e.g. sickle cell trait)
 - patients that have received a blood transfusion within the past 3 weeks
 - patients that have received cancer chemotherapy within the past 3 weeks
- In cases of rapidly evolving type 1 diabetes the increase of HbA1c values might be delayed compared to the acute increase in glucose concentrations. In these conditions diabetes mellitus must be diagnosed based on plasma glucose concentration and/or the typical clinical symptoms.
- Diluted samples cannot be used with Afinion HbA1c.
- Coagulated or hemolysed samples cannot be used with Afinion HbA1c.
- If the sample has a hemoglobin value below 6,0 g/dL or above 20,0 g/dL, no test result will be reported and an information code will be displayed.

QUALITY CONTROL

Quality control testing should be done to confirm that your Afinion System is working properly and providing reliable results. Only when controls are used routinely and the values are within acceptable ranges can accurate results be assured for patient samples.

 It is recommended to keep a permanent record of all quality control results. The Afinion Analyzer automatically stores the control results in a separate log. Consult the Afinion Analyzer User Manual.

Choosing control material

 Afinion HbA1c Control from Abbott is recommended for routine quality control testing. Consult the Afinion HbA1c Control Package Insert.

If controls from another supplier are used, the precision must be determined and acceptable ranges for the Afinion System must be established. Lyophilised cell free lysate or hemolysed whole blood control materials cannot be used with Afinion HbA1c.

Frequency of control testing

Controls should be analysed:

- anytime an unexpected test result is obtained.
- with each shipment of Afinion HbA1c Test Kits.
- with each new lot of Afinion HbA1c Test Kits.
- when training new operators in correct use of the Afinion HbA1c and the Afinion Analyzer.
- in compliance with national or local regulations.

Verifying the control results

The measured value should be within the acceptable limits stated for the control material. Consult the Afinion HbA1c Control Package Insert.

If the result obtained for the Afinion HbA1c Control is outside the acceptable limits, make sure that:

- patient samples are not analysed until control results are within acceptable limits.
- the control vial has not passed its expiry date.
- the control vial has not been used for more than 60 days.
- the control vial and Afinion HbA1c Test Cartridge have been stored according to recommendations.
- there is no evidence of bacterial or fungal contamination of the control vial.

Correct any procedural error and retest the control material.

If no procedural errors are detected:

- examine the laboratory's quality control record to investigate the frequency of control failures.
- ensure that there is no trend in out-of-range quality control results.
- retest the control material using a new control vial.
- Patient results must be declared invalid when controls do not perform as expected. Contact your local supplier for advice before analysing patient samples.

TROUBLESHOOTING

To ensure that correct HbA1c results are reported, the Afinion Analyzer performs optical, electronic and mechanical controls of the capillary, the test cartridge and all individual processing steps during the course of each analysis. When problems are detected by the built-in failsafe mechanisms, the Analyzer terminates the test and displays an information code.

The table below contains Afinion HbA1c specific information codes. Consult the Afinion Analyzer User Manual for information codes not listed in this table.

Code #	Cause
103	The hemoglobin concentration is below 6,0 g/dL
104	The hemoglobin concentration is above 20,0 g/dL
105	The HbA1c value is below measuring range
106	The HbA1c value is above measuring range

PERFORMANCE CHARACTERISTICS

The performance data presented in this section are representative data from internal and external studies. Results obtained in individual laboratories may vary.

Method comparison

Fingerstick and venous whole blood samples from 120 patients (4.6-11.4 % HbA1c), 38-42 at each of three sites, were analysed with the Alere Afinion AS100 Analyzer using three different Afinion HbA1c lots. The venous samples were sent to a laboratory for duplicate analysis with an HPLC method. The results are shown in Table 1 and Table 2.

Table 1: Method comparison. Afinion HbA1c vs. a laboratory HPLC method. Weighted Deming regression slope and intercept for 120 fingerstick samples (3 sites, 3 lots).

	Estimate	95 % lower bound	95 % upper bound
slope	1.00	0.97	1.03
Intercept (% HbA1c)	0.00	-0.21	0.20

Table 2: Method comparison. Afinion HbA1c vs. a laboratory HPLC method. Weighted Deming regression slope and intercept for 120 venous whole blood samples (3 sites, 3 lots).

	Estimate	95 % lower bound	95 % upper bound
slope	0.99	0.96	1.02
Intercept (% HbA1c)	0.05	-0.16	0.27

Precision

Within device (total) precision of Afinion HbA1c used with the Alere Afinion AS100 Analyzer was determined according to CLSI Protocol EP5-A2. Afinion HbA1c Control C I, Control C II, and four EDTA whole blood samples were analysed for 10 days, 2 runs per day and 4 replicates per run. 3 lots of Afinion HbA1c Test were used. Precision data are summarised in Table 3 (% HbA1c) and Table 4 (HbA1c, mmol/mol).

Table 3: Within device (total) precision. N (Number of replicates per lot)=80. CV=Coefficient of Variation.

Sample	Lot 1		Lot 2		Lot 3	
	% HbA1c	CV (%)	% HbA1c	CV (%)	% HbA1c	CV (%)
Low	5.2	1.1	5.1	1.1	5.0	1.3
Threshold	6.5	1.2	6.4	1.1	6.4	0.9
Medium	8.4	1.1	8.4	1.1	8.4	1.0
High	11.8	1.0	11.8	0.9	11.8	1.2
Control C I	6.3	1.0	6.1	1.0	6.2	0.9
Control C II	8.3	0.9	8.5	0.8	8.4	0.8

Table 4: Within device (total) precision. N (Number of replicates per lot)=80. CV=Coefficient of Variation.

Sample	Lot 1		Lot 2		Lot 3	
	HbA1c (mmol/mol)	CV (%)	HbA1c (mmol/mol)	CV (%)	HbA1c (mmol/mol)	CV (%)
Low	33	2.0	32	1.6	32	1.9
Threshold	47	1.7	46	1.6	46	1.4
Medium	69	1.5	68	1.6	68	1.2
High	106	1.2	106	1.1	105	1.4
Control C I	46	1.4	43	1.5	44	1.5
Control C II	68	1.4	69	1.0	68	1.0

Performance testing with the AFINION™ 2 Analyzer

The performance of Afinion HbA1c used with the Afinion 2 Analyzer has been demonstrated to be equivalent to the performance obtained with the Alere Afinion AS100 Analyzer.

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1. Bunn F et al. The Biosynthesis of Human Hemoglobin A1c. Slow glycosylation of hemoglobin in vivo. *J Clin Invest* 1976; 57:1652-1659.
2. International Expert Committee, International Expert Committee Report on the Role of the A1c Assay in the Diagnosis of Diabetes. *Diabetes Care* 2009; 32:1327-1334.
3. The Diabetes Control and Complications Trial Research Group, The Effect of Intensive Treatment of Diabetes on the Development and Progression of Long-Term Complications in Insulin-Dependent Diabetes Mellitus. *N Engl J Med* 1993; 329:977-986.
4. Sacks DB et al., Guidelines and Recommendations for Laboratory Analysis in the Diagnosis and Management of Diabetes Mellitus. *Clin Chem* 2002; 48:436-472.
5. Stratton IM et al., Association of glycemia with macrovascular and microvascular complications of type 2 diabetes: prospective observational study (UKPDS 35). *BMJ* 2000; 321:405-412.
6. Jeppsson, JO et al., Approved IFCC Reference Method for the Measurement of HbA1c in Human Blood, *Clin Chem Lab Med* 2002; 40(1):78-89.
7. The American Diabetes Association, European Association for the Study of Diabetes, International Federation of Clinical Chemistry and Laboratory Medicine. Consensus Statement on the Worldwide Standardization of the Hemoglobin A1C Measurement. *Diabetes Care* 2007; 30:2399-2400.
8. Weykamp C et al., The IFCC Reference Measurement System for HbA1c: a 6-year progress report. *Clin Chem* 2008; 54: 240-248.
9. Hoelzel W et al., IFCC reference system for measurement of haemoglobin A1c in human blood and the national standardization schemes in the United States, Japan, and Sweden: a method-comparison study. *Clin Chem* 2004; 50: 166-174.
10. Nathan DM et al., A1c-Derived Average Glucose Study Group. Translating the A1C assay into estimated average glucose values. *Diabetes Care* 2008; 31:1473-1478

Abbott Diagnostics Technologies AS

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P.O. Box 6863 Rodeløkka

NO-0504 Oslo, Norway

www.globalpointofcare.abbott

ISO 13485 certified company.

© 2021 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners.