## RESULTADOS

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- Stone NJ, Robinson JG, Lichtenstein AH, et al. 2013 ACC/AHA guideline on the treatment of blood cholesterol to reduce atherosclerotic cardiovascular risk in adults: A Report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines. *Circulation* 2014; 129:S1-S45.
- Bachorik PS, Ross JW, for the National Cholesterol Education Program Working Group on Lipoprotein Measurement. National Cholesterol Education Program recommendations for measurement of low-density lipoprotein cholesterol: executive summary. *Clin Chem* 1995;41:1414-20.
- Goff DC, Lloyd-Jones DM, Bennett G, et al. 2013 ACC/AHA guideline on the assessment of cardiovascular risk: A Report of the American College of Cardiology/American Heart Association Task Force on Practice Guidelines. *Circulation* 2014; 129:S49-S73.
- Tietz NW, ed. *Fundamentals of Clinical Chemistry*. Philadelphia, Pa.: WB Saunders Co., 1987.
- Siedel J, Hagele EQ, Ziegenhorn J, et al. Reagent for the enzymatic determination of serum total cholesterol with improved lipolytic efficiency. *Clin Chem* 1983;29:1075-80.
- Warnick GR, Benderson J, Albers JG. Dextran sulfate-Mg2+ precipitation procedure for quantitation of high density lipoprotein cholesterol. *Selected Methods for Clinical Chemistry* 1983;10:91-9.
- Allain CC, Poon LS, Chan CS, et al. Enzymatic determination of total serum cholesterol. *Clin Chem* 1974;20:470-5.
- Roeschlau P, Bernt E, Gruber W. Enzymatische bestimmung des gesamtcholesterins im serum. *Z Klin Chem Klin Biochem* 1974;12:226.
- Fossati P, Prencipe L. Serum triglycerides determined colorimetrically with an enzyme that produces hydrogen peroxide. *Clin Chem* 1982;28:2077-80.
- Jessen RH, Dass CJ, Eckfeldt JH. Do enzymatic analyses of serum triglycerides really need blanking for free glycerol? *Clin Chem* 1990;36:1372-5.
- Tietz NW, ed. *Clinical Guide to Laboratory Tests*. 2nd ed. Philadelphia, Pa.: WB Saunders Co., 1990.
- Blumenfeld TA, Hertelendy WG, Ford SH. Simultaneously obtained skin- puncture serum, skin-puncture plasma, and venous serum compared, and effects of warming the skin before puncture. *Clin Chem* 1977;23:1705-10.
- CV Risk Calculator. http://professional.heart.org/professional/GuidelinesStatements/PreventionGuidelines/UCM\_457698.Prevention-Guidelines.jsp, accessed 5/26/17.
- Jacobson TA, Ito MK, Maki KC, et al. National Lipid Association recommendations for patient-centered management of dyslipidemia: Part 1 – executive summary. *J Clin Lipidol* 2014;8:473-88.
- Castelli WP, Abbott RD, McNamara PM. Summary estimates of cholesterol used to predict coronary heart disease. *Circulation* 1983;67:730-4.
- Kinosian B, Glick H, Garland G. Cholesterol and coronary heart disease: predicting risks by levels and ratios. *Ann Intern Med* 1994;121:641-7.
- American Diabetes Association. Standards of medical care in diabetes—2017, 2. Classification and Diagnosis of Diabetes. *Diabetes Care* 2017;40(Suppl.1):S11-24.

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