



## Safety Data Sheet (SDS)

acc. to OSHA HazCom Standard 29 CFR 1910.1200(g) and Canada WHMIS 2015

Issue date 10/29/2018

Reviewed on 10/29/2018

### \* 1 Product and Company Identification

- **Product Identifier**
- **Trade Name:** Wash Buffer
- **Index number:** MW2
- **Relevant identified uses of the substance or mixture and uses advised against:**  
No further relevant information available.
- **Application of the substance / the mixture:** Kit Component
- **Details of the Supplier of the Safety Data Sheet:**
- **Manufacturer/Supplier:**  
Company, Manufacturer:  
TECHLAB, Inc.  
2001 Kraft Drive  
Blacksburg, VA 24060, USA

Distributed by:

Alere North America, LLC

30 South Keller Road

Orlando, Florida 32810, US

- **Emergency telephone number:** +(207)730-5750

### 2 Hazard(s) Identification

- **Classification of the substance or mixture:**



Skin Sens. 1B H317 May cause an allergic skin reaction.

- **Hazard description:** Refer to the information found in this Safety Data Sheet.
- **Information concerning particular hazards for human and environment:**  
The product has to be labeled due to the calculation procedure of international guidelines.
- **Classification system:**  
The classification was made according to the latest editions of international substances lists, and expanded upon from company and literature data.

- **Label elements:**

- **GHS label elements**

The product is classified and labeled according to the Globally Harmonized System (GHS).

- **Hazard pictograms:**



GHS07

- **Signal word:** Warning

- **Hazard-determining components of labeling:**

reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1)  
Sodium Dodecyl Sulphate

(Contd. on page 2)

## Safety Data Sheet (SDS)

acc. to OSHA HazCom Standard 29 CFR 1910.1200(g) and Canada WHMIS 2015

Issue date 10/29/2018

Reviewed on 10/29/2018

### Trade Name: Wash Buffer

#### • Hazard statements:

H317 May cause an allergic skin reaction.

#### • Precautionary statements:

P261 Avoid breathing dust/fume/gas/mist/vapors/spray.

P272 Contaminated work clothing must not be allowed out of the workplace.

P280 Wear protective gloves.

P302+P352 If on skin: Wash with plenty of water.

P362+P364 Take off contaminated clothing and wash it before reuse.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

P321 Specific treatment (see supplementary first aid instructions on this Safety Data Sheet).

P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

#### • Unknown acute toxicity:

This value refers to knowledge of known, established toxicological or ecotoxicological values.

0 % of the mixture consists of component(s) of unknown toxicity.

#### • Hazard(s) not otherwise classified (HNOC): None known

#### • Results of PBT and vPvB assessment

• **PBT:** Not applicable.







• **vPvB:** Not applicable.

### 3 Composition/Information on ingredients:

#### • Chemical characterization: Mixtures

• **Description:** Mixture of substances listed below with non-hazardous additions.

#### • Dangerous Components:

CAS: 9002-93-1	Polyethylene glycol octylphenol ether  	0.5%
CAS: 151-21-3 EINECS: 205-788-1 RTECS: WT 1050000	Sodium Dodecyl Sulphate  	0.1%
CAS: 55965-84-9 Index number: 613-167-00-5	reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1)  	0.05%

#### • SVHC

9002-93-1 Polyethylene glycol octylphenol ether

• **Additional information:** For the wording of the listed hazard phrases refer to section 16.

### 4 First-Aid Measures

#### • Description of first aid measures

#### • General information:

Symptoms of poisoning may occur after exposure to dust, fumes or particulates; seek medical attention if feeling unwell.

(Contd. on page 3)

## Safety Data Sheet (SDS)

acc. to OSHA HazCom Standard 29 CFR 1910.1200(g) and Canada WHMIS 2015

Issue date 10/29/2018

Reviewed on 10/29/2018

**Trade Name: Wash Buffer**

- **After inhalation:**  
Remove to fresh air and keep at rest in a position comfortable for breathing. Obtain medical attention if breathing difficulty persists.
- **After skin contact:**  
Remove contaminated clothes and shoes.  
Immediately wash with water and soap and rinse thoroughly.  
If skin irritation occurs, consult a doctor.
- **After eye contact:**  
Rinse opened eye for several minutes under running water.  
If easy to do so, remove contact lenses if worn.  
Seek immediate medical advice.
- **After swallowing:**  
Call a doctor immediately.  
Seek immediate medical advice.  
Do not induce vomiting without medical advice.
- **Information for doctor**
- **Most important symptoms and effects, both acute and delayed:**  
Prolonged exposure may cause an allergic reaction.
- **Indication of any immediate medical attention and special treatment needed:**  
No further relevant information available.

### 5 Fire-Fighting Measures

- **Extinguishing media**
- **Suitable extinguishing agents:**  
CO<sub>2</sub>, extinguishing powder or water spray. Fight larger fires with water spray.  
Use fire fighting measures that suit the environment.
- **For safety reasons unsuitable extinguishing agents:** Water with full jet
- **Special hazards arising from the substance or mixture:** No further relevant information available.
- **Advice for firefighters**
- **Special protective equipment for firefighters:** Mouth respiratory protective device.

### 6 Accidental Release Measures:

- **Personal precautions, protective equipment and emergency procedures:**  
Do not breathe vapor.  
Avoid contact with skin, eyes and clothing.
- **Environmental precautions:** Do not allow to enter sewers/surface or ground water.
- **Methods and material for containment and cleaning up:**  
Dispose contaminated material as waste according to section 13.  
Ensure adequate ventilation.  
Dispose of the collected material according to regulations.
- **Reference to other sections:**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protection equipment.  
See Section 13 for disposal information.

### 7 Handling and Storage

- **Handling**
- **Precautions for safe handling:**  
Avoid contact with skin, eyes and clothing

(Contd. on page 4)

## Safety Data Sheet (SDS)

acc. to OSHA HazCom Standard 29 CFR 1910.1200(g) and Canada WHMIS 2015

Issue date 10/29/2018

Reviewed on 10/29/2018

### Trade Name: Wash Buffer

- Wear protective equipment.
- Avoid creating and breathing dust/fume/gas/mist/vapors/spray.
- Ensure good ventilation/exhaustion at the workplace.
- **Information about protection against explosions and fires:** No special measures required.
- **Conditions for safe storage, including any incompatibilities**
- **Storage**
- **Requirements to be met by storerooms and receptacles:**
  - Protect from heat and direct sunlight.
  - Store in a cool, dry place.
  - Store in the original container.
- **Information about storage in one common storage facility:** Not required.
- **Further information about storage conditions:** Keep receptacle tightly sealed.
- **Specific end use(s):** No further relevant information available.

### 8 Exposure Controls/Personal Protection

- **Additional information about design of technical systems:** No further data; see section 7.
- **Control parameters:**

#### Components with occupational exposure limits:

##### 56-81-5 Glycerol

PEL (USA)	Long-term value: 15* 5** mg/m <sup>3</sup> mist; *total dust **respirable fraction
TLV (USA)	TLV withdrawn-insufficient data human occup. exp.
EL (Canada)	Long-term value: 10* 3** mg/m <sup>3</sup> *mist; **mist, respirable
EV (Canada)	Long-term value: 10 mg/m <sup>3</sup>

- **Additional information:** The lists that were valid during the creation of this SDS were used as basis.
- **Exposure controls:**
- **Personal protective equipment**
- **General protective and hygienic measures:**
  - Keep away from foodstuffs, beverages and feed.
  - Immediately remove all soiled and contaminated clothing and wash before reuse.
  - Wash hands before breaks and at the end of work.
- **Breathing equipment:**
  - Have self-contained breathing apparatus, positive pressure, MSHA/NIOSH (approved or equivalent) available in case of spillage or equipment failure.
- **Protection of hands:**



Protective gloves

- **Material of gloves:** Nitrile rubber, NBR
- **Penetration time of glove material:** Not applicable.

(Contd. on page 5)

## Safety Data Sheet (SDS)

acc. to OSHA HazCom Standard 29 CFR 1910.1200(g) and Canada WHMIS 2015

Issue date 10/29/2018

Reviewed on 10/29/2018

**Trade Name: Wash Buffer**

• **Eye protection:**



Safety goggles

• **Body protection:**



Protective work clothing

### 9 Physical and chemical properties

• **Information on basic physical and chemical properties**

• **General Information**

• **Appearance:**

**Form:** Liquid

**Color:** Clear

• **pH-value:** Not applicable.

• **Change in condition**

**Melting point/Melting range:** Not determined.

**Boiling point/Boiling range:** Not determined.

• **Flash point:** Not applicable.

• **Flammability (solid, gaseous):** Not determined.

• **Ignition temperature:** Not applicable

• **Decomposition temperature:** Not determined.

• **Auto igniting:** Product is not self-igniting.

• **Danger of explosion:** Product does not present an explosion hazard.

• **Explosion limits:**

**Lower:** Not determined.

**Upper:** Not determined.

• **Vapor pressure:** Not applicable.

• **Density:**

**Relative density:** Not determined.

**Vapor density:** Not applicable.

**Evaporation rate:** Not applicable.

• **Solubility in / Miscibility with:**

**Water:** Miscible

• **Partition coefficient (n-octanol/water):** Not determined.

• **Viscosity:**

**Dynamic:** Not applicable.

**Kinematic:** Not applicable.

(Contd. on page 6)

## Safety Data Sheet (SDS)

acc. to OSHA HazCom Standard 29 CFR 1910.1200(g) and Canada WHMIS 2015

Issue date 10/29/2018

Reviewed on 10/29/2018

**Trade Name:** Wash Buffer

• **Solvent content:**

• **Solids content:**

0.0 %

• **Other information:**

No further relevant information available.

### 10 Stability and reactivity

- **Reactivity:** Stable under normal conditions.
- **Chemical stability:** Stable under normal conditions.
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- **Possibility of hazardous reactions:** No dangerous reactions known.
- **Conditions to avoid:**  
Incompatibilities, flames, ignition sources.  
High temperatures.  
Keep out of direct sunlight.
- **Incompatible materials:** Strong acids, strong bases, strong oxidizing agents and strong reducing agents.
- **Hazardous decomposition products:** Oxides

### 11 Toxicological information

• **Information on toxicological effects:**

• **Acute toxicity:**

• **LD/LC50 values that are relevant for classification:**

**56-81-5 Glycerol**

Oral	LD50	12,600 mg/kg (Rat)
Dermal	LD50	>10,000 mg/kg (Rabbit)
Inhalative	LC50/96 hours	>50,000 mg/l (Trout)

**9002-93-1 Polyethylene glycol octylphenol ether**

Oral	LD50	1,800 mg/kg (Rat)
Dermal	LD50	8,000 mg/kg (Rat)
Inhalative	LC50/4 h	8.9 mg/l (Pimephales)

**151-21-3 Sodium Dodecyl Sulphate**

Oral	LD50	1,288 mg/kg (Rat)
Inhalative	LC50/4 h	>3,900 mg/l (Rat)
THE LC50 VALUES HEREUNDER ARE ESTIMATED ON THE BASIS OF A 4-HOUR EXPOSURE.		

• **Primary irritant effect:**

• **On the skin:** May be irritating.

• **On the eye:** May be irritating.

• **Sensitization:** Sensitization possible through skin contact.

• **Additional toxicological information:**

The product shows the following dangers according to internally approved calculation methods for preparations:

Irritant

• **OSHA-Ca (Occupational Safety & Health Administration):**

None of the ingredients are listed.

(Contd. on page 7)

## Safety Data Sheet (SDS)

acc. to OSHA HazCom Standard 29 CFR 1910.1200(g) and Canada WHMIS 2015

Issue date 10/29/2018

Reviewed on 10/29/2018

**Trade Name: Wash Buffer**

### 12 Ecological information

• **Toxicity:**

• **Aquatic toxicity:**

**9002-93-1 Polyethylene glycol octylphenol ether**

EC50 26 mg/l (Daphnia)

• **Persistence and degradability:** No further relevant information available.

• **Behavior in environmental systems:**

• **Bioaccumulative potential:** No further relevant information available.

• **Mobility in soil:** No further relevant information available.

• **Additional ecological information:**

• **General notes:**

Do not allow undiluted product or product that has not been neutralized to reach ground water, water course or sewage system.

• **Results of PBT and vPvB assessment:**

• **PBT:** Not applicable.

• **vPvB:** Not applicable.

• **Other adverse effects:** No further relevant information available.

### 13 Disposal considerations

• **Waste treatment methods**

• **Recommendation:**

Observe all federal, state and local environmental regulations when disposing of this material.

• **Uncleaned packaging**

• **Recommendation:** Dispose of as unused product.

### 14 Transport information

• **UN-Number:**

• **ADR/ADN, ADN, IMDG, IATA**

Non-Regulated Material

• **UN proper shipping name:**

• **ADR/ADN, ADN, IMDG, IATA**

Non-Regulated Material

• **Transport hazard class(es):**

• **ADR/ADN, ADN, IMDG, IATA**

• **Class:**

Non-Regulated Material

• **Packing group:**

• **ADR/ADN, IMDG, IATA**

Non-Regulated Material

• **Environmental hazards:**

Not applicable.

• **Special precautions for user:**

Not applicable.

• **Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code:**

Not applicable.

• **UN "Model Regulation":**

Non-Regulated Material

(Contd. on page 8)



## Safety Data Sheet (SDS)

acc. to OSHA HazCom Standard 29 CFR 1910.1200(g) and Canada WHMIS 2015

Issue date 10/29/2018

Reviewed on 10/29/2018

**Trade Name: Wash Buffer**

### 15 Regulatory information

- **Safety, health and environmental regulations/legislation specific for the substance or mixture:**
- **SARA (Superfund Amendments and Reauthorization):**

- **Section 355 (extremely hazardous substances):**

None of the ingredients are listed.

- **Section 313 (Specific toxic chemical listings):**

None of the ingredients are listed.

- **TSCA (Toxic Substances Control Act):**

7647-14-5 Sodium Chloride

56-81-5 Glycerol

5704-04-1 N-(tri(hydroxymethyl)methyl)glycine

9002-93-1 Polyethylene glycol octylphenol ether

151-21-3 Sodium Dodecyl Sulphate

- **TSCA new (21st Century Act) (Substances not listed)**

55965-84-9 reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1)

- **Chemicals known to cause cancer:**

None of the ingredients are listed.

- **Chemicals known to cause reproductive toxicity for females:**

None of the ingredients are listed.

- **Chemicals known to cause reproductive toxicity for males:**

None of the ingredients are listed.

- **Chemicals known to cause developmental toxicity:**

None of the ingredients are listed.

- **New Jersey Right-to-Know List:**

56-81-5 Glycerol

- **New Jersey Special Hazardous Substance List:**

None of the ingredients are listed.

- **Pennsylvania Right-to-Know List:**

56-81-5 Glycerol

- **Pennsylvania Special Hazardous Substance List:**

None of the ingredients are listed.

- **Carcinogenic categories:**

- **EPA (Environmental Protection Agency):**

None of the ingredients are listed.

- **TLV (Threshold Limit Value established by ACGIH):**

None of the ingredients are listed.

- **NIOSH-Ca (National Institute for Occupational Safety and Health):**

None of the ingredients are listed.



## Safety Data Sheet (SDS)

acc. to OSHA HazCom Standard 29 CFR 1910.1200(g) and Canada WHMIS 2015

Issue date 10/29/2018

Reviewed on 10/29/2018

**Trade Name:** Wash Buffer

• **Canadian substance listings:**

• **Canadian Domestic Substances List (DSL):**

All ingredients are listed.

• **Canadian Ingredient Disclosure list (limit 0.1%):**

None of the ingredients are listed.

• **Canadian Ingredient Disclosure list (limit 1%):**

None of the ingredients are listed.

• **GHS label elements**

The product is classified and labeled according to the Globally Harmonized System (GHS).

• **Hazard pictograms:**



GHS07

• **Signal word:** Warning

• **Hazard-determining components of labeling:**

reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1)

Sodium Dodecyl Sulphate

• **Hazard statements:**

H317 May cause an allergic skin reaction.

• **Precautionary statements:**

P261 Avoid breathing dust/fume/gas/mist/vapors/spray.

P272 Contaminated work clothing must not be allowed out of the workplace.

P280 Wear protective gloves.

P302+P352 If on skin: Wash with plenty of water.

P362+P364 Take off contaminated clothing and wash it before reuse.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

P321 Specific treatment (see supplementary first aid instructions on this Safety Data Sheet).

P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

• **National regulations:**

None of the ingredients are listed.

• **Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

\* **16 Other information:**

The information and recommendations in this safety data sheet are, to the best of our knowledge, accurate as of the date of issue. Nothing herein shall be deemed to create warranty, expressed or implied, and shall not establish a legally valid contractual relationship. It is the responsibility of the user to determine applicability of this information and the suitability of the material or product for any particular purpose.

• **Relevant phrases:**

H301 Toxic if swallowed.

H302 Harmful if swallowed.

H311 Toxic in contact with skin.

H314 Causes severe skin burns and eye damage.

H315 Causes skin irritation.

(Contd. on page 10)

## Safety Data Sheet (SDS)

acc. to OSHA HazCom Standard 29 CFR 1910.1200(g) and Canada WHMIS 2015

Issue date 10/29/2018

Reviewed on 10/29/2018

### Trade Name: Wash Buffer

- H317 May cause an allergic skin reaction.
- H319 Causes serious eye irritation.
- H331 Toxic if inhaled.
- H335 May cause respiratory irritation.
- H400 Very toxic to aquatic life.
- H410 Very toxic to aquatic life with long lasting effects.
- R20/22 Harmful by inhalation and if swallowed.
- R22 Harmful if swallowed.
- R23 Toxic by inhalation.
- R23/24/25 Toxic by inhalation, in contact with skin and if swallowed.
- R34 Causes burns.
- R36 Irritating to eyes.
- R36/37/38 Irritating to eyes, respiratory system and skin.
- R43 May cause sensitisation by skin contact.
- R50/53 Very toxic to aquatic organisms, may cause long-term adverse effects in the aquatic environment.

### Abbreviations and acronyms:

- ADR: The European Agreement concerning the International Carriage of Dangerous Goods by Road
- ADN: The European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways
- IMDG: International Maritime Code for Dangerous Goods
- IATA: International Air Transport Association
- ACGIH: American Conference of Governmental Industrial Hygienists
- EINECS: European Inventory of Existing Commercial Chemical Substances
- ELINCS: European List of Notified Chemical Substances
- CAS: Chemical Abstracts Service (division of the American Chemical Society)
- LC50: Lethal concentration, 50 percent
- LD50: Lethal dose, 50 percent
- PBT: Persistent, Bioaccumulative and Toxic
- SVHC: Substances of Very High Concern
- vPvB: very Persistent and very Bioaccumulative
- NIOSH: National Institute for Occupational Safety and Health
- OSHA: Occupational Safety & Health Administration
- Acute Tox. 4: Acute toxicity – Category 4
- Acute Tox. 3: Acute toxicity – Category 3
- Skin Corr. 1B: Skin corrosion/irritation – Category 1B
- Skin Irrit. 2: Skin corrosion/irritation – Category 2
- Eye Irrit. 2A: Serious eye damage/eye irritation – Category 2A
- Skin Sens. 1: Skin sensitisation – Category 1
- Skin Sens. 1B: Skin sensitisation – Category 1B
- STOT SE 3: Specific target organ toxicity (single exposure) – Category 3
- Aquatic Acute 1: Hazardous to the aquatic environment - acute aquatic hazard – Category 1
- Aquatic Chronic 1: Hazardous to the aquatic environment - long-term aquatic hazard – Category 1

### \* Data compared to the previous version altered.

© 2018 TECHLAB, Inc. All Rights Reserved

The TECHLAB Logo and TECHLAB are trademarks of TECHLAB, Inc., under license.

The Alere Logo is a trademark of the Abbott group of companies.