

	<p style="text-align: center;">External Release Notes Abbott Rapid Diagnostics (ARDx) Informatics</p>		<p>Template No. D10036037 Rev AA</p>
<p>Document Type:</p>	<p>FORM</p>	<p>Governing Document:</p>	<p>D10046628</p>

1. Introduction

Welcome to RALS™ version 8.4.1 build 23.43.0342. This document details the fix found in this product version. For information regarding the installation, features, functionality, and use of this product, please refer to the Resources section.

2. System Requirements

- For existing customers remaining on an ARDx Informatics Provided Image (R70 platform) the following minimal requirements should be available to upgrade to this version of RALS:
 - C Drive must have at least 10 Gb of free space.
 - D Drive must have three times the size of the backup file.
- For customers implementing in a Customer Provided Environment the following minimal requirement should be available:
 - See [Specifications for Customer Provided Environments](#)

3. Supported Environments

- Supported Browsers:
 - Microsoft Edge
 - Google Chrome
 - Mozilla Firefox

4. Required Third Party Components

- No additional third-party components required. All components are installed on the RALS server as part of the product installation.

5. Bug Fixes

- The Abbott i-STAT® result values of ACT WBT received from the device are converted to ACT-C and ACT-K as appropriate based on the result test type for transfer to LIS/EMR system.

6. Resources

www.rals.com

7. Contact Information

If you have any questions regarding this information, contact Abbott Rapid Diagnostics Informatics at:

rals.support@abbott.com or 877-627-7257